KIMTECH™

Kimtech™ A5 Sterile Overboots



Features **CLEAN-DON™ technology** for simpler aseptic donning

Made from **low-lint, breathable SMS material** with a cloth-like feel

Sterile overboots with **anti-slip soles**

Kimtech™ A5 Sterile Overboots improve worker safety and comfort, while maintaining cleanliness, providing seamless protection when and where it counts. The aseptic donning technique saves time while reducing the risk of contamination, and the overboots are suitable for use in EU GMP ISO 5 Grade A sterile cleanrooms.

The overboots feature an anti-slip vinyl sole and two extra-long ties for a secure fit. Simple, fast and secure donning is ensured through the boots' elastic opening and their inside-out folding. The sterile cleanroom overboots are made from high performance SMS material which has a cloth-like feel. This low-lint fabric features high-strength bound and triple-stitched seams and adheres to Helmke Drum test Category I. The gamma-irradiated overboots also feature a high bacterial filtration efficiency.

The overboots provide a comfortable and secure fit and are delivered pair vacuum-packed for added sterility assurance. This efficient packaging is also designed to reduce space so more products can fit into the gowning area storage.

KIMTECH™

Kimtech™ A5 Sterile Overboots

Light Duty



Long Boots



Key Features

- Featuring innovative CLEAN-DON™ technology for fast, simple and compliant aseptic donning
- Contamination protection is ensured through high quality, sterility assured material with high-strength bound and triple-stitched seams
- Pair vacuum-packed

Assured Compliance

- PPE Cat I according to Regulation (EU) 2016/425
- Certificate of Conformance and Certificate of Sterility by lot available online

Quality Standards

- Suitable for EU GMP ISO 5 Class A sterile cleanroom environments
- Manufactured in accordance with Quality System ISO 9001
- Sterility Assurance Level 10^{-6}
- Particle cleanliness level Helmke Drum test Category I
- Sterilisation in accordance to ISO 11737-1, -2

Applications

- EU GMP ISO 5 Grade A
- Aseptic Processing
- Parental Drug Manufacturing
- Biotechnology
- Pharmaceutical Compounding
- Ophthalmic Product Manufacturing

CE

Kimtech™ A5 Sterile Overboots			
Code	**Description**	**Size**	**Total/Case**
12922	Light Duty	S (29.2cm)	100 pairs
88808	Light Duty	Universal (34.9cm)	100 pairs
31683	Long Boots	S (29.2cm)	100 pairs
31696	Long Boots	Universal (34.9cm)	100 pairs