Product Specification Sheet

VWR PURE WHITE ROUND BOROSILICATE GLASSES OF HYDROLYTIC CLASS I

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release.The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Ø (mm)	Pk
631-1576	10	1.000
631-1577	12	1.000
631-1578	13	1.000
631-1579	15	1.000
631-1580	18	1.000
631-1581	20	1.000
631-1582	22	1.000
631-1583	24	1.000
631-1584	25	1.000
631-1585	30	1.000

Thickness: 0.13 - 0.16 mm

Material: made of chemically resistant borosilicate glass D 263™ M of the first hydrolytic class

Other Characteristics:

- high spectral transmission
- virtually colourless appearance
- excellent optical quality
- low flatness deviations
- good chemical resistance
- refractive index fine-tuned to microscopes
- low autofluorescence

Product Specification Sheet

- protective coating facilitates automated separation of individual coverglasses in coverslippers
- complies precisely with the set thickness tolerances
- persistency at changing climatic conditions
- low alkaline strength
- No bubbles, no flaws
- Pre-washed

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **July 2014**