epredia Enhancing precision cancer diagnostics

Technical Data Sheet

Product Identification	**Microscope slides, Non Adhesion, for in-vitro-diagnostic use**

Article code 1 – 18											
No. Area	**1**	**2**	**3-8**	**9**	**10**	**11**	**12**	**13**	**14-16**	**17**	**18**
Art. code digits	**A**	**X**	**0000XX**	**X**	**X**	**E**	**0**	**X**	**XXX**	**X**	**X**
Description	Product Type-Microscope Slides	Edge Finish	Dimension / Colors	Tab	Thickness	Glass Type-soda lime	Box Type-50 pcs cardboard box	Cello Definition	Private Label Epredia Label	Carton Size	Packaging

Usage	ready to use, single use only
Raw material	extra white glass
Refractive Index	n (λ = 546 nm): 1,53 ± 0,02
Abbe value	n/a

Dimensions	**ISO 8037-1:** length: 76,0 mm (+0 / -1 mm) width: 26,0 mm (+0 / -1 mm)
Thickness	**1,0 mm** = 0,95 - 1,05 mm **1,1 mm** = 1,0 - 1,20 mm

Intended Use & Quality Assurance	
Standards and directive:	See Declaration of Conformity
Quality standards:	1. cleanliness 2. dimensions 3. damages 4. adhesion / stuck 5. wettability
Storage and usage instructions:	See Instruction for use

This document is valid without an original signature. Please contact New Erie Scientific LLC (Quality Management), a subsidiary of Epredia, in any case of questions.

Epredia
New Erie Scientific LLC
20 Post Road | Portsmouth, NH 03801, United States
Tel: +1 (603) 431-8410