ThermoFisher SCIENTIFIC

EP Scientific Products, LLC
Product Specification Sheet

Catalog Number:	**130-04A/PF**
Description:	EP Scientific Products, LLC Certified 4 oz amber short straight side wide mouth jar with closure.
Cap:	EP#-00357
Color	White
Size	58-400
Solid or Septa Style	Solid
Material of Construction	Polypropylene
Cap Liner:	CLT1009-58C-MJ
Color	Natural
Material(s) of Construction	PTFE-faced Silicone
Bonded or Unbonded	Bonded
Container:	EP#-01164
Color	Amber
Style	4 oz Wide-Mouth Short Straight Side Bottle
Dimensions D x H x ID	2.25 x 2.75in (57.15 x 69.85mm)*
Thread Finish	58-405
Material of Construction	Type III Soda-Lime Glass
Processing Level:	Certified
Certification Parameters:	Endotoxin <0.06 EU/mL
Manufacturing Address:	520 N. Main Street. Miami, OK 74354

*Dimensions are nominal and subject to change without notice.