Triethanolamine
Baker Analyzed

Material No: 7078.9025
Effective from: 2014/08/13

Product Specification

Test	Specification
Assay	99.0–103.0 %
Appearance of solution	passes test
Density (g/ml) at 25°C	1.117–1.125
Heavy Metals (as Pb)	<= 10 ppm
Identification	passes test
Iron (Fe)	<= 5 ppm
Sulfated Ash	<= 0.005 %
Water (H2O)	<= 0.2 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000