Hospidex
Safety solutions



Hospidex

sharpsafe

PRODUCT DATA SHEET

General data	
Item code	SHA-51106610 (without flap) - SHA-51105610 (with flap)
GTIN code	5056271508436 (without flap) - 5056271508429 (with flap)
Product Name	Sharpsafe® 5th GEN. 7 liters (flap in option)
Description	Sharps container 7 liters componed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste

Characteristics data	
Nominal volume	7 l
Total volume	6,76 l
Usefull volume	5,54 l
Maximum weight	1,35 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	175 mm
Wide	175 mm
Height	380 mm
Opening	55 x 55 mm
Fill line height	310 mm
Weight	401 g (410 g avec clapet)
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 100% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	62%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availibility Accessories	Stands, trolley brackets (see brochure and compatibility tables)



Logistic data	
Carton quantity	50
Carton weight	23,5 Kg
Carton dimensions (LxWxH)	810 x 345 x 680 mm
Carton volume	190,0262 dm3
Palet quantity	450 units
Pallet weight	227 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm

Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes

Manufacture data	
Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group