Identification of the product

Product number	3806001
Product type	Shell vials
GMDN - Code	none
GMDN - Definition	none
HS-Code	70179000
UDI-DI (smallest sales unit)	04250317332592
UDI-DI (Multipack)	

Manufacturer / Supplier Identification

Trademark owner	Paul Marienfeld GmbH & Co.KG Am Wöllerspfad 4, 97922 Lauda-Königshofen, Germany
Country of origin	Germany

Specification & Compliance

Norm	
Dimensions	35x12 mm
Material	soda lime glass
Type of neck or rim	straight rim
Wall thickness	thickness of wall 0.8 - 1.0 mm
Packing unit per box	144 pcs./carton
Brand	MARIENFELD SUPERIOR

Packaging

Smallest sales unit	144 pieces
Multipack	