
CORNING

Corning Life Sciences
836 North Street
Building 300 Suite 3401
Tewksbury MA 01876

1/7/2015

Dear Valued Customer,

This letter is in response to your request for the dimensions on Falcon® Catalog # 351006, 50mmx9mm Bacteriological Petri Dish with Tight-Fit Lid, Polystyrene.



All dimensions in millimeters (mm)

ITEM	CAT. NO.	DESCRIPTION	DISH AND LID KEY DIMENSIONS											
			A	B	C	D	E	F	G	H	J	K	L	M
1	353005	50x9 mm Bacteriological Dish. Non-treated. 20/pack, 500/case	50.82	48.93	48.32	50.27	55.06	8.26	7.17	1.09	52.92	51.4	4.26	0.99

Please feel free to contact me at DLQAHelp@corning.com with any questions.

Sincerely,

Patricia Clarke

Patricia Clarke
Regulatory Compliance & Quality Systems Manager
Corning Life Sciences - Discovery Labware
DLQAhelp@corning.com