Specification

Buffer Solution pH 12.454 277089

Physical Description:	Liquid
Product Code:	277089
Product Name:	Buffer Solution pH 12.454
Short Description:	usage: These buffer solutions are useful for the calibration of high precision laboratory pH meters
Headline Comment:	* standardized, according to DIN19266
Specifications:	pH (20°C): 12.627 ± 0.01 pH (10°C): 12.810 ± 0.01 pH (25°C): 12.454 ± 0.01 pH (30°C): 12.289 ± 0.01
WGK:	nwg
Storage:	RT
CS:	38221900

AppliChem GmbH
Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • **info.de@itwreagents.com** • **www.itwreagents.com**
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340