Technical Data Sheet

Article number	**40567002**
Description	**Pasteur pipettes**

Material	Soda Lime Glass
Length	230 mm
Form	With elongated tip
Wall	0.55 mm
Body	approx. 90 x 7 mm Ø
Outer Ø of capillary tip	approx. 1.5 mm
Volume	approx. 1.9 ml

ISO 7712
Uniform, exact execution, single use

Tel.: +49 9779/808-0
Fax: +49 9779/808-88
Stettener Str. 22-24
97647 Sondheim/Rhön · Germany
E-Mail: info@assistent.eu
Internet: http://www.assistent.eu

Glaswarenfabrik Karl Hecht GmbH & Co KG | Amtsgericht Schweinfurt HRA 8063 | USt.IdNr./VAT-No.: DE132198535
Persönlich haftende Gesellschafterin: Hecht GmbH | AG Schweinfurt HRB1160
Geschäftsführer: Karl-Heinz Hecht, Thomas Diepholtz | Einzelprokura: Tonja Müller