VITLAB
Competence in Labware

Sample preparation, Wash-bottles

VITsafe™ safety wash bottles, narrow-mouth

ASK A QUESTION

KNOWLEDGE: VITSAFE™



Bottles made from PE-LD (or PP for Acetone and MEK), spray tube from PP.

With coloured VENT-CAP which supports easy visual identification of the media.

More safety due to the durable safety imprint in accordance to Directive (EC) No. 1272/2008 (GHS), as well as with all important information:
- Material name in German, English, French and Spanish
- Chemical formula, CAS No., hazard pictogram, signal word
- Risk phrases (H phrases), safety phrases (P phrases), as well as NFPA Code

With precise spray jet and optimised medium backflow. Practically no leakage or dripping with the VENT-CAP, the design of which prevents almost all static overpressure.

Types / Sizes

Imprint	Colour VENT-CAP	Volume ml	Material	Thread	Ø mm	Height mm	PU	Cat. No.
Acetone	red	250	PP	GL 25	65	135	12	1431829
Acetone	red	500	PP	GL 25	74	180	12	1432829
Acetone	red	1000	PP	GL 32	92	215	12	1433829
Acetonitrile	red	500	PE-LD	GL 25	74	180	6	1332969
Distilled Water	blue	250	PE-LD	GL 25	65	135	12	1331819
Distilled Water	blue	500	PE-LD	GL 25	74	180	12	1332819
Distilled Water	red	1000	PE-LD	GL 32	92	221	12	1333819
Acetic acid	red	500	PE-LD	GL 25	74	180	6	1332979

Imprint	Colour VENT-CAP	Volume ml	Material	Thread	Ø mm	Height mm	PU	Cat. No.
Ethanol	orange	250	PE-LD	GL 25	65	135	12	1331869
Ethanol	orange	500	PE-LD	GL 25	74	180	12	1332869
Ethanol	red	1000	PE-LD	GL 32	92	221	12	1333869
Ethyl acetate	red	250	PE-LD	GL 25	65	135	12	1331859
Ethyl acetate	red	500	PE-LD	GL 25	74	180	12	1332859
Ethyl acetate	red	1000	PE-LD	GL 32	92	221	12	1333859
Isopropanol	yellow	250	PE-LD	GL 25	65	135	12	1331849
Isopropanol	yellow	500	PE-LD	GL 25	74	180	12	1332849
Isopropanol	red	1000	PE-LD	GL 32	92	221	12	1333849
Methanol	green	250	PE-LD	GL 25	65	135	12	1331839
Methanol	green	500	PE-LD	GL 25	74	180	12	1332839
Methanol	red	1000	PE-LD	GL 32	92	221	12	1333839
Methylene chloride	transparent	500	PE-LD	GL 25	74	180	6	1332879
Methyl ethyl ketone (MEK)	red	500	PP	GL 25	74	180	6	1432989
Heptane	red	500	PE-LD	GL 25	74	180	6	1332899
i-Hexane	red	500	PE-LD	GL 25	74	180	6	1332909
N,N-Dimethylformamide	red	500	PE-LD	GL 25	74	180	6	1332889
Pentane	red	500	PE-LD	GL 25	74	180	6	1433959
Tetrahydrofuran (THF)	red	500	PE-LD	GL 25	74	180	6	1332939
Toluene	red	500	PE-LD	GL 25	74	180	6	1332949
Xylene	red	500	PE-LD	GL 25	74	180	6	1332959

Height without spray tube

Other variations available upon request.

You might also be interested in



VENT-CAP wash bottle caps, PP



Wash bottles with imprint, PE-LD/PP



VITsafe™ safety wash bottles, wide-mouth

© 2023 VITLAB GmbH www.vitlab.com. All rights reserved.