phone 1.847.816.5070
fax 1.847.816.5072
info@heathrowscientific.com

Heathrow Scientific®

620 Lakeview Parkway
Vernon Hills, IL 60061 USA
www.heathrowscientific.com

TECHNICAL DATA SHEET

Description of Quality Standard: Thermoformed PVC with non-toxic white pigment which is inert or minimally reactive with most common laboratory reagents. Inspected in the molding facility for cracks, incomplete, distorted shape, damage or defects which would compromise proper function. Sealed in air-tights bags. Manufactured in a clean room with sterile versions certified as RNase/DNase and pyrogen free.
Product is not manufactured or derived from and does not come in contact with material of animal origin.
Certified phthalates, latex and asbestos free.
Sterile versions are sterilized by Gamma radiation.

Product name:	Solution Basin, 50 mL, Non-Sterile, PVC
Part number:	120690
Color:	White
Material of construction	Thermoformed PVC