Glycerol
Baker Analyzed

Material No: 7044.9025
Effective from: 2013/06/28

Product Specification

Test	Specification
Assay	>= 99.5 %
Color	passes test
Density (g/ml) at 25°C	>= 1.255
Fatty Acid Esters	<= 0.05 %
Neutrality	passes test
Residue after Ignition	<= 0.005 %
Silver Reducing Substances	passes test
Substances Darkened by H2SO4	passes test
Sulphate (SO4)	<= 0.001 %
Trace Impurities (in ppm)	
Chloride (Cl)	<= 5
Heavy Metals (as Pb)	<= 2

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000