Product Specification

Material	85041.001
Material description	Buffer pH 4.00 ± 0.02 (20°C)
Grade	AVS TITRINORM Reagent traceable to SRM from NIST

Additional information

Characteristics	Specifications
pH (20°C) (tolerance ± 0.02)	3.98 - 4.02
NIST traceable	Confirmed

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001 Leuven; Belgium

Additional information

The pH of this buffer solution is traceable to and verified against primary Standard Reference Materials (SRM) from National Institute of Standards and Technology (NIST): SRM 188 and SRM 186 I + II.

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.