Grade 1289

Item number 207

Characteristics:

Qualitative filter paper, hardened, plain

Applications:

Medium fast filtering

Physical properties	Units	Typical Value
Basis Weight	g/m^2	**84**
Filtration Speed	s / 10 ml	**20**
Burst Strength, wet	kPa	**≥30**
Typical Retention (*)	µm	**8-12**

18.12.2017 Rev02

(*)Approximate value, not part of the product specification.

All indicated data to be understood as typical average values.

This product is manufactured in conformance with the Quality Management standard DIN EN ISO 9001.

Contact office:

Ahlstrom-Munksjö Germany GmbH
Niederschlag 1
D-09471 Bärenstein
Germany

T: +49 (0) 37347 / 83-0
F: +49 (0) 37347 / 83-64
Email: filtrak@ahlstrom-munksjo.com

DISCLAIMER: The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as Ahlstrom-Munksjö deems appropriate. Refer to www.ahlstrom-munksjo.com for contact information.

© Ahlstrom-Munksjö 2017