Product Specification Sheet

VWR FOIL ROLLS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release.The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Thickness	Pk
113-1576	Foil rolls, 102 mm×153 m	0.015 mm +/- 5%	24
113-1577	Foil rolls, 153 mm×153 m	0.015 mm +/- 5%	24

Material: Aluminium
Temperature resistance : 200°C to -660°C
Standards : 21 CR 178.3910 Ultraviolet light absorbance requirement

Other Characteristics:

- Small size allows for use in limited spaces and easy transport
- Designed for an easy use and minimal waste.
- Custom-sized pieces of foil can easily be created for covering small surfaces such as beakers, flasks, tubes and petri dishes.
- Ideal for autoclaving, freezing, incubating, protecting lab ware, lining, sealing and general storage.
- Designed exclusively for the Mini Bin™ Foil Dispenser.
- Formed of aluminum foil and rolled on cardboard roller.

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **June 2014**