Product Specification Sheet

VWR FORMALIN NEUTRALISER PADS, Q PATH® SAFETY PAD

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release.The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Dimensions	Pk
720-2251	Q Path® Safety Pad	203×254 mm	25
720-2252	Q Path® Safety Pad	280×356 mm	25
720-2253	Q Path® Safety Pad	406×485 mm	25

Other Characteristics:

- Absorbs the odours (neutralises formalin).
- Ensures good contrast for biopsies.
- Protects knives and cutting boards against premature wear.
- Different sizes available

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **August 2014**