
CATALOG NUMBER	CLOSURE SIZE	MATERIALS				
		CLOSURE	GASKET	TUBING	PORT CAP	VENT CAP
2162-0531	53B (1/4" PORT)	PP (WHITE)	TPE	NALGENE 50 PLATINUM CURED SILICONE	TPE	TPE

THREAD FINISH (mm).....:	53	
THREAD (T.P.I.).....:	4	
NOMINAL WEIGHT (ASSEMBLED gm).....:	62	
	INCH	MM
A (NOMINAL HEIGHT).....:	2.7	69
D (DIA.).....:	2.59	65.5
BARBED PORT SIZE (INCH)...:	1/4	
VENT PORT SIZE (INCH).....:	1/4	

NOTE:

53B CLOSURES ARE SUPPLIED WITH A 4" SHORT TUBE AND A 18" LONG TUBE ATTACHED.



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV.
PB	GRP	PSS 00-2162-0531	0413
CK'D KOW	4/14/13		

DESCRIPTION:

53B FILLING/VENTING CLOSURE