TECHNICAL DATA SHEET

29076935

CONT-DAIRY---PP---CLEA---3,0--

3 OZ VIAL

Single Use Flip-Top Vial for Dairy, Food Sampling



Applications

Intended Fields of Use

- Field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials

Options

- Can be supplied with label (barcode…)
- Standard vial is colorless – colors upon request

Characteristics

Characteristics	Typical Values
Internal Volume	77 ml
Filling Level Line	60 ml
Vial Overall Height	85.5 mm
Top Body Diameter	36.3 mm
Bottom Body Diameter	35.1 mm
Vial Weight	11.1 grams

Delivery

Qty Per Box	400 units
Qty Per Pallet	12000units

Benefits

Flip-Top Vial Design

- User-friendly hinged vials concept, one-hand opening and closing
- Patented “in mold closing” process
- Translucent in color for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use
- Allows the use of gloves in cold environment

Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

(Rev. 06/2024)

Please contact us for official regulatory submission documents.