n-Hexane 99%

ChromAR HPLC

Material No: 6713-25
Effective from: 2013/03/01

Product Specification

Test	Specification
Assay (sum of 5 isomers, total hexanes, plus methylcyclopentane) (by GC, corrected for water)	>= 99 %
Assay (as n-hexane (by GC, corrected for water)	>= 95.0 %
Water (KF)	<= 0.01 %
Water Soluble Titrable Acid (meq/g)	<= 0.0005
Titrable Base (meq/g)	<= 0.0003
Residue after Evaporation	<= 0.0005 %
Ultraviolet Absorbance (1.00-cm cell vs. water):	
210 nm	<= 0.2
220 nm	<= 0.05
230 nm	<= 0.02
240 nm	<= 0.009

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000