Type in Product Names, Product Numbers, or CAS Numbers to see suggestions.

Specification Sheet

Product Name	Mayer’s Mucicarmine Stain Solution for microscopy
Product Number	41325
Product Brand	SIAL
CAS Number	51395-97-2
Molecular Weight	

TEST	SPECIFICATION
APPEARANCE (COLOR)	Purple to Dark Purple
APPEARANCE (FORM)	Liquid
APPEARANCE (TURBIDITY)	Clear
SUITABILITY ASSAY	FUNGAL SPORES AND HYPHAE APPEAR PALE TO DEEP PURPLE (OBSERVED UNDER MICROSCOPE BY USING OIL IMMERSION LENSE)
EXPIRY PERIOD	36 MONTHS

Feedback

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Reproduction of any materials from the site is strictly forbidden without permission.
Site Use Terms | **Privacy Policy** | **General Terms and Conditions of Sale** | **Copyright Consent** | **Cookies Settings**