Spécifications du produit

Produit	22715.293
Description du produit	Chloroforme
Grade / qualité	SPECTRONORM Stabilized with about 1% of ethanol
Numéro CAS	67-66-3
Formule moléculaire	$CHCl_3$

Information additionnelle

Analyses	Spécifications
Titre (GC)	≥ 99.0 %
Acidité	≤ 0.0002 meq/g
Alcalinité	≤ 0.0002 meq/g
Coloration	≤ 10 APHA
Résidu à l'évaporation	≤ 2 ppm
Eau	≤ 100 ppm
Absorbance (245 nm)	≤ 0.82
Absorbance (250 nm)	≤ 0.30
Absorbance (255 nm)	≤ 0.22
Absorbance (260 nm)	≤ 0.07
Absorbance (à partir de 270 nm)	≤ 0.01
Transmittance (245 nm)	≥ 15 %
Transmittance (250 nm)	≥ 50 %
Transmittance (255 nm)	≥ 60 %
Transmittance (260 nm)	≥ 85 %
Transmittance (à partir de 270 nm)	≥ 98 %
Fluorescence (en quinine) (254 nm)	≤ 1.0 ppb
Fluorescence (en quinine) (365 nm)	≤ 1.0 ppb

Signature

Nous certifions que ce lot est conforme aux spécifications listées ci-dessus.

Ce document a été produit électroniquement et est validé sans signature.

VWR International

Information additionnelle

Répond aux exigences de la "Spectroscopie à usage spécifique" selon ACS

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.

VWR International LLC, Radnor Corporate Center, Building One, Suite 200, 100 Matsonford Road, Radnor, PA 19087, USA
VWR International bv, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium