Lindenstraße 43 - 97892 Kreuzwertheim
Tel. + 49 (0) 9342 9283-0
Fax +49 (0) 9342 39860
Internet: www.amarell.de
E-Mail: **thermometer@amarell.de**

Technical Data Sheet

Article number:	E904010
Type:	Digital thermometer type “Speed Lab”
Measurement range:	-50+300°C
Resolution:	0,1°C
Accuracy:	+/-0,5°C from -10 up to +100°C +/-1,5°C from +100,1 up to 199,9°C +/-2,5°C from 200 up to +300°C
Battery:	Round cell, CR2032
Battery life:	c. 5700 h
Dimensions:	c. 155x26mm
Operation temp.:	0 to +50°C
Probe:	Semi conductor probe made of stainless steel
Response time:	T_{90} = 2,5 seconds
Probe length:	75x2,5x1,6mm
Special functions:	Max/min function, hold function, water resistant, °C ↔ °F

January 2022