Whirl-Pak® Specification Sheet
B01323WA

Size: 7 ½" W x 15" L (19 x 38 cm).

Capacity: 69 oz. (2041 ml).

Film: Blend of low-linear and low-density virgin polyethylene, 3.0 mil (.076 mm) thick, FDA approved for food contact.

Closing Tab: Aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality system certified to **ISO9001.**

Packaging: 500 bags/box in a liner bag twist-tied shut.

Condition of Use: Not to be used at temperatures above 180º F (82º C).

Phone 1.800.558.9595 | Fax 1.800.372.1236 | Website whirl-pak.com | E-mail info@whirl-pak.com
901 Janesville Avenue, Fort Atkinson, WI 53538-0901