
3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Safranin O – certified by the Biological Stain Commission

Product Number:	**S8884**
CAS Number:	477-73-6
MDL:	MFCD00011759
Formula:	C20H19ClN4
Formula Weight:	350.84 g/mol



TEST	Specification
Appearance (Color)	Green to Very Dark Green
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Red to Very Dark Red
Solubility (Form)	Solution
1 mg/mL H2O	
Biological Stain Commission Certificatio	Pass
Commission certified number: AcS-46	
Dye Content	≥ 80 %

Specification: PRD.1.ZQ5.10000064751

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.