
# Technical Product Data

Technical Data Sheet | Labcon Product Number 3181-345-008-9




Dimensions are inches. Do not scale





Image is not actual size. Colors may differ.

## 50 mL SuperClear® Centrifuge Tubes, 50 per Bag, Sterile

| | | | |
|---|---|---|---|
| **Fluid Measure** | mL | **Freezable** | Yes |
| **Fluid Volume** | 50 | **Freeze Temp** | -80ºC |
| **Height (mm)** | 115 | **Autoclavable** | 121ºC 15min |
| **Width (mm)** | 29 | **Sterile** | Yes |
| **Length (mm)** | | **Sterile Processing** | Certified as sterile with a sterility assurance level (SAL) of 10-6 |
| **Pack Qty** | 50 | **Centrifugation** | Yes |
| **Case Qty** | 500 | **Centrifugation Speed (RCF)** | 20000 |
| **Package Type** | Bag | **Shelf Life (Months)** | 36 |
| **Material** | Polypropylene | **Certification** | |
| **Cap Material** | Polypropylene | **Regulatory** | |

**Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.**

**Additional Testing:** ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

IATA 95kPa certified to conform with USDOT HMR, 49 C.F.R., Parts 171-180 and IATA DGR UN3373 PI650 mailing tube standards and packaging instructions for diagnostic and infectious specimens.

This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is **© Copyright 03-May-2021** . All Rights Reserved. Technical Data Sheet Version 1/2020