Crucible

Iron
Area of use: Only on open burner flame.

Article no.	crucible Ø mm	height mm	thickness mm	
5620	40	40	1,5	55 g
5621	45	40	1,5	63 g
5622	48	36	1,5	78 g
5623	50	45	2,0	94 g
5624	60	50	1,5	106 g



Cover	cover Ø mm	
5630	40	16 g
5631	45	18 g
5632	48	22 g
5633	50	22 g
5634	60	28 g

© 2015

Carl Friedrich Usbeck KG
Industriestr. 12-14
42477 Radevormwald
Germany

Phone: +49 (02195) 91 18-0
Telefax: +49 (02195) 91 18-40

https://www.usbeck.eu