BAKERBOND spe, Octadecyl (C18)

BAKERBOND spe

Disposable Extraction Columns; Glass

Material No: 7334-04
Effective from: 2015/07/07

Product Specification

Test	Specification
The following results are based on the bulk bonded phase:	
Appearance	passes test
Capacity Factor: Butylbenzene	For Information Only
Capacity Factor: Ethylbenzene	For Information Only
Capacity Factor: Propylbenzene	For Information Only
Carbon (C)(%)	For Information Only
Desisopropyl Atrazine (10%), ml	110 - 150
Surface Coverage (ug/m2), Based on %C	320 - 350
The following results are based on the silica gel raw material:	
Average Particle Diameter, um (APD) - via Malvern (100 mm lens)	47 - 60
The following results are based on the solid phase extraction columns:	
Extraction Residue	<= 0.2 %

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
1000 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000