45280 Lithium 6, ^{6}Li, plasma standard solution, Specpure®, ^{6}Li 100µg/ml

Product Number:	45280
MDL number:	MFCD00084621
Molecular formula:	Matrix: 2% HNO_3

Product Specification

Comment:	This is a non-radioactive, "stable" isotope.
Concentration:	100 µg/ml 6 Lithium
Concentration:	100µg/mL ±2%
Matrix:	2% HNO3

Date of Print:	September 6, 2022
Version:	1

Product Specifications are subject to amendment and may change over time.