BENZA-CLEAN®

Cleansing pads for external use

HOS-10-5201HD



Description
- Fan folded pads sized 20 x 13 cm saturated with 0,113% benzalkonium chloride - contains no alcohol.
- Hermetically sealed in a laminated aluminium foil pouch
- For single and external use only

Application
- Cleansing pads for professional, hospital & private use.

Commercial Data

• Item ID	HOS-10-5201HD
• Type	BENZA-CLEAN® - Cleansing pads for external use
• Basic UDI-DI	542501868BENZACLEAN0000KZ

Technical Data

• Material	Fan folded pad saturated with 0,113% benzalkonium chloride
• Validity	5 years
• Marking	CE-marked, Class 1, non-sterile Lot number and expiration date printed on each pouch, box and case
• Material Safety Data Sheet	Available upon request
• Certificate of Analysis	Available with every batch produced upon request
• Quality Management	Hospidex NV / Hospidex France SAS are ISO EN 13485:2016 certified

Packaging

	Box - 100 towelettes	Case - 20 boxes	Pallet - 60 cases
• Dimensions	18 x 10 x 6 cm	41 x 22 x 34 cm	120 x 80 x 215 cm
• Weight	0,360 kg	8 kg	495 kg