Scheda tecnica

Technical data sheet

SUPPORTO ASSORBENTE "SPONGY"/ *"SPONGY" ADSORBENT HOLDER*

cod./cat. PBIB88733.A, PBIB24762.A, PBIB88734.A





CARATTERISTICHE TECNICHE

Il supporto assorbente "Spongy" é realizzato in poliuretano espanso, color antracite.

Antiurto, protegge le provette (contenitori primari) da cadute accidentali limitando le possibilità di rottura.

Funziona anche da assorbente in caso di rotture accidentali quando inserito nei contenitori secondari dedicati durante il trasporto.

"Spongy" adsorbent holder is made of expanded polyurethane, anthracite colour.

Shock-Proof, protects the test tubes (primary containers) from accidental drops reducing the possibility of break.

It works as an adsobent sponge in case of accidental breaks when used inside the dedicated secondary containers.

COME ORDINARE / HOW TO ORDER

Cod./Cat.	**Tipo di contenitore secondario/Type of Secondary Container**	**n. di provette da 16/13 mm/N. of test tubes 16/13 mm**	**Dimensioni h mm / Dimensions h mm**	**Pezzi per cf/ Pcs per case**
PBIB24762.A	"Biocarrier"	80/100	310x170x50	10
PBIB88733.A	"Bio-Protection"	80/100	300x125x30	10
PBIB88734.A	"Leak-Proof"	150/180	360x210x30	1