

REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
C	REVISED DIMENSIONS TO MATCH NEW MOLD	JEU	08/05/19	201906041106214





NECK DETAIL

SCALE: 2 X 1

NOTES:
1. 1 1/2 TURNS OF THREAD.
2. ALL DIMENSIONS APPLY AFTER FLAMING.
3. MINIMUM WALL THICKNESS: .015.
4. WEIGHT: 6.2 ±0.5 GRAMS.
5. CAPACITY: 37.5 ±2.5 C ML.

ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
3	144	20mm CAP | MATERIAL: WHITE POLYPROPYLENE
2	144	20mm TIP | MATERIAL: NATURAL LINER LOW DENSITY POLYETHYLENE
1	144	30mL BOTTLE | MATRERIAL: WHITE LOW DENSITY POLYETHYLENE

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
DIMENSIONS IN INCHES UNLESS SPECIFIED
C .XX ± .01"
.XXX ± .015"
.XXXX ± .0005"
FRACTIONAL ± 1/64"

THIRD ANGLE PROJECTION

DURAN WHEATON KIMBLE

TITLE: 30ML BTL, DPG, LDPE WHT, ECDT

ITEM NUMBER:	DRAWING NUMBER:	DRAWN BY:	DRAWN DATE:
211624	CSD-211624	BL	07/07/05
SAP NUMBER: 211624	OLD MATERIAL NUMBER: N/A	SHEET SIZE: B	SHEET NUMBER: 1 OF 1 / SHEET SCALE: 1:1