BAKERBOND spe, PAH Aqua

BAKERBOND spe

Disposable Extraction Columns

Material No: 7490-08
Effective from: 2015/10/23

Product Specification

Test	Specification
Amino (NH2):	
Carbon (C)(%)	For Information Only
Extraction Residue	<= 0.2 %
Nitrogen (N)(%)	For Information Only
Particle Size Distribution: APD, microns	47 - 61
Specific Surface Area, m2/g	For Information Only
Surface Coverage (ug/m2), Based on %C	For Information Only
Surface Coverage (ug/m2), Based on %N	For Information Only
Oktadecyl (C18):	
Appearance	passes test
Average Particle Diameter, um (APD) - via Malvern (100 mm lens)	47 - 60
Capacity Factor: Butylbenzene	For Information Only
Capacity Factor: Ethylbenzene	For Information Only
Capacity Factor: Propylbenzene	For Information Only
Carbon (C)(%)	For Information Only
Desisopropyl Atrazine (10%), ml	110 - 150
Extraction Residue	<= 0.2 %
Surface Coverage (ug/m2), Based on %C	320 - 350

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
1500 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000