
:\GCA\W217004

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF WHEATON INDUSTRIES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF WHEATON INDUSTRIES INC
ALL RIGHTS RESERVED ©
ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF WHEATON'S SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

THIRD ANGLE PROJECTION

PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION
C 1.	24	16oz BOTTLE
C 2.	24	70-400 CAP W/ LINER
B 3.	24	50077397 LABEL (NOT SHOWN)



BOTTLE MATERIAL: TYPE III CLEAR SODA-LIME GLASS
CAP MATERIAL: WHITE POLYPROPYLENE
B D LINER MATERIAL: 1.02 [.040] POLYETHYLENE FOAM WITH 0.13 [.005] PFTE-FACED

REV	E.C.N. #	DATE
D	2014012714492591	02/27/14 TP
	REVISED CAP VIEW MAX/MIN WAS ±	9/18/13 TP
	REVISED CAP VIEW	03/14/13 TP
C	2013011411102976	01/16/13 TP
B	2011101211223320	11/03/11 TP
A	20070131111503	02/06/07 BL

4. DIMENSIONS NOT TOLERANCED ARE FOR REFERENCE ONLY. D
3. SIGHT "H" DOWN "E".
2. DEPRESSED THREADS AT SEAMS.
1. DIMENSIONS ARE IN MILLIMETERS AND [INCHES]. B

----NOTES----

TITLE: 16oz BTL, W/M AC RD, TYPE III CLR, PTFE LN B
CAPACITY: 16oz (517.5cc ±7.4cc OVERFLOW)
FINISH: G.P.I. 70-400-405 B
SCALE: FULL UNLESS NOTED, DO NOT SCALE THIS DRAWING.

DATE DRAWN: 02/06/07
DRAWN BY: B. LLOYD
CHECKED BY: TRJ 2/7/07
PROJ #:
APPROX. WGT.: 9.473oz (268.6gr)
AUTHORIZED: CJM 2/7/07

WHEATON®
1501 North 10th Street | Millville, NJ 08332-2038 | USA

DWG. NO: B- W217004

UNCONTROLLED COPY
DEC 17 2015


K:\GCA\W216947

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF WHEATON INDUSTRIES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF WHEATON INDUSTRIES INC. ALL RIGHTS RESERVED.© ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF WHEATON'S SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

THIRD ANGLE PROJECTION




ISOMETRIC VIEW
SCALE 1 :2



B C BOTTLE MATERIAL: TYPE III AMBER SODA-LIME GLASS
CAP MATERIAL: WHITE POLYPROPYLENE
LINER MATERIAL: .040 [1.02] THICK POLYETHYLENE FOAM WITH .005 [0.13] PTFE

PRODUCT COMPONENTS		
ITEM	QTY	DESCRIPTION
B 1.	24	4oz BOTTLE
B 2.	24	38-400 CAP W/LINER
3.	24	LABEL (NOT SHOWN)



FINISH DETAIL
SCALE 2 : 1

5.	APPROX WEIGHTS: BOTTLE 118.7gr. | CAP 2.8±1gr. EXCLUDING LINER. D
4.	BOTTOM VIEW MAY VARY DUE TO DIFFERENT MANUFACTRING LOTS. D
3.	DEPRESSED THREADS AT SEAMS. SIGHT "H" AT "E".
2.	MILLIMETER DIMENSIONS AND DIMENSIONS NOT TOLERANCED ARE FOR REFERENCE ONLY.
1.	DIMENSIONS ARE IN INCHES AND [MILLIMETERS]. B
	----NOTES----

TITLE: 4oz BTL, WM PKR, TYPE III AMB, PP/PTFE LN B
DATE DRAWN: 05/05/04 | PROJ#:
DRAWN BY: B. LLOYD | APPROX WGT: SEE NOTE #5
CAPACITY: 137-143cc OVERFOLW
CHECKED BY: TRJ 05/12/04 | AUTHORIZED: CJM 05/12/04
FINISH: G.P.I. 38-400-405 C/R B | SIZE: 120cc PRACTICAL FILL
SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.

WHEATON®
1501 North 10th Street, Millville, NJ 08332-2038

SIZE: B
DWG.NO: W216947

REV	ECN #	DATE
D	2015032608311340	04/27/15 TP
C	201401270935495	01/28/14 TP
B	2011101211223320	11/02/11 TP
A	20040407081917	05/05/04 BL

UNCONTROLLED COPY
DEC 17 2015