DWK Life Sciences Limited

Unit 4, Riverside 2, Campbell Road
Stoke-on-Trent, ST4 4RJ, UK

Registered in England and Wales No.7723574

Tel: +44 (0)1782 444406
Fax: +44 (0)1782 940436

www.dwkltd.com
enquiries.uk@dwk.com

PRODUCT SPECIFICATION

Specification	**Product Code**	360251
	Description	250ml Container, No Label, Metal Screw Cap
	Material Base	Polystyrene
	Material Cap	Metal
	Colour	Clear
	Sterility	Sterilised by ETO
	Shelf Life	5 Years from Manufacture

Dimension	**Height (mm)**	121.2
	Diameter (mm)	63.4
	Packaging Format	25/Pack 50/Case

Compliance	**Food safe (EU/FDA)**	EU 10/2011 / FDA CFR 21 177
	BSE / TSE free	YES - No tallow from animal sources are used in the manufacture of the finished product.
	REACH Compliant	Compliant according to 1907/2006
	LATEX Free	YES - No Materials containing latex natural rubber are used in the manufacturing, handling, and packaging process.
	CE Marking	Meets the essential requirements of European directive 98/79 EC
	Leak Tested	In Accordance with BS EN 14254 Annexe D and 95KPa compliant in accordance with ASTM D4991 Standard

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: R. Adams
Product Manager

Date: 15th March 2021