Methanol
ULTRA RESI-ANALYZED
For Purge and Trap Analysis

Material No.: 9077-02
Revision No.: 0

Product Specification

Test	Specification
Assay (CH_3OH) (by GC, corrected for water)	≥ 99.9 %
Residue after Evaporation	≤ 1.0 ppm
Titrable Acid (µeq/g)	≤ 0.3
Titrable Base (µeq/g)	≤ 0.10
Water (by KF, coulometric)	≤ 0.08 %
Volatile Organic Trace Analysis - Below EPA 8260B CRQL	Conforms

For Laboratory,Research,or Manufacturing Use
Performance Tested for Use in EPA Methods
500 Series for Drinking Water
600 Series for Wastewater
846 for Solid Waste

Packaging Site: Phillipsburg Mfg Ctr & DC

Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Mansford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700