NB600-381

Orders: orders@novusbio.com
Support: technical@novusbio.com
Web: www.novusbio.com

Protocols, Publications, Related Products, Reviews and more:
www.novusbio.com/NB600-381

NB600-381

V5 Epitope Tag Antibody

Unit Size
0.1 ml

Concentration
1.0 mg/ml

Storage
Store at 4C. Do not freeze.

Applications
WB, ELISA, ICC/IF, IP

Species
Non-species specific

Host
Rabbit

Clonality
Polyclonal

Purity
Immunogen affinity purified

Preservative
0.09% Sodium Azide

Updated 3/9/2013 2.0

Immunogen
Rabbits were immunized with a synthetic peptide representing amino acid residues 95 to 108 of RNA polymerase alpha subunit of simion virus 5 conjugated to KLH.

Recommended Dilutions
ELISA, Immunocytochemistry/Immunofluorescence 1:100- 1:400, Immunoprecipitation 1 - 4 ug/mg lysate, Western Blot 1:1000-1:10000


Buffer
Phosphate Buffered Saline (PBS)

Uses
This antibody is useful for ELISA, Western blot, and immunocytochemistry.


References (more available at www.novusbio.com/NB600-381)
Hance MW, Dole K, Gopal U et al Secreted Hsp90 is a novel regulator of the epithelial to mesenchymal transition (EMT) in prostate cancer J Biol Chem 2012 Sep 18 [PMID: 22989880] (WB

Images

Western Blot: V5 Epitope Tag Antibody [NB600-381] - NP-40 whole cell lysate (10 ug) from wild-type HEK293 cells (lane 1) or HEK293 cells expressing V5 tagged Nicastrin (lane 2). NB600-381 was used at 0.2 mcg/ml (1:5000).



Western Blot: V5 Epitope Tag Antibody [NB600-381] - 200, 100, or 50 ng of E. coli whole cell lysate expressing a multi-tag fusion protein. Antibody used at 0.04 ug/ml (1:25,000).



Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. ***For more information on our guarantee, please visit www.novusbio.com/guarantee.***

Novus USA
8100 Southpark Way, A-8
Littleton, CO 80120
p) 888-506-6887
p) 303-730-1950
f) 303-730-1966
novus@novusbio.com

Novus Canada
461 North Service Road West,
Unit B37
Oakville, ON L6M 2V5
p) 855-668-8722
f) 905-827-6402
canada@novusbio.com

Novus Europe
12 Cambridge Science Park
Cambridge, UK CB4 0FQ
UK: p) +44-1223-426-001, f) +44-871-971-1635
DE: p) +49-6922-22340-60, f) +49-0800-58926-79
IT: p) +39-02-4032-6786, f) +39-02-4032-6340
FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com