



## KLINIPATH Microscope Slides

| | |
|---|---|
| Description: | Clear imprinted with KLINIPATH LOGO ; Clear with a "★" at each of the two lower corners; pre-cleaned, clear; no marks, prints or contamination, etc.<br>Color: Standard colors<br>Color Frosted Area: 19.5 x 25 (±0.5)mm; |
| Purpose: | Used for histological pathology examination |
| Materials: | Standard Glass; thinckness 1.0 (-0.05,+0.1)mm |
| Features: | The smooth edges provide the safety of the operation |
| | Color frosted |
| | Suitable for marking with conventional labeling system, pencil, mark pens and thermal transfer printers |
| | Silane coated, to provide permanent adhesion which can bind the tissue sections to the glass |
| Temperature: | Normal laboratory temperature and humidity |
| Uses: | Recommended for routine H&E stains,IHC,ISH,frozen sections |
| Shelf life： | 18 Months |

### KLINIPATH Microscope Slides

| EU Part No. | Dimension | | | | Frosted or Painted | Paint Color | Corner Angle | Edges | Packaging | QTY/CASE |
|---|---|---|---|---|---|---|---|---|---|---|
| | A (mm) | B (mm) | C (mm) -0.05/+0.1 | D (mm) (+/-0.5) | | | | | | |
| KP-SIL-3056 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | WHITE | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / | 2500 |
| KP-SIL-3057 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | BLUE | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / | 2500 |
| KP-SIL-3058 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | YELLOW | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / | 2500 |
| KP-SIL-3059 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | PINK | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / | 2500 |
| KP-SIL-3060 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | GREEN | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / | 2500 |
| KP-SIL-3061 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | ORANGE | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / Case | 2500 |
| KP-SIL-3062 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | Lilac | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / Case | 2500 |
| PR-S-001 | 75.7~76.0 | 25.0~25.3 | 1 | N/A | Painted | WHITE | 45° (small corner) | Smooth | 50pcs / Box, 50 Boxes / Case | 2500 |