
CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS	
		CARBOY	CLOSURE
002214	5	HDPE	HDPE

SIZE CODE. . . . . :	-0050
SIZE (NOM VOL L/GAL):	20/5
NOM BRIM CAP (L). . . . :	20.2
NECK FINISH (mm). . . . :	70
THREAD (T.P.I.). . . . . :	4
NOM. WEIGHT W/ CLOSURE(gms). . . . . :	2300
A (NECK I.D.). . . . . :	2.32 ±.05
B (OVERALL HEIGHT). . :	15.6 ±.4
D (WIDTH). . . . . :	12.56 ±.25
E (DEPTH). . . . . :	9.08 ±.20
W.T.(NOM WALL THK). . :	0.10



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Therno Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV.
MCC	PVC 2/15/91	PSS 00-2214-0050	—
CK'D PVC			

DESCRIPTION:

DOT-34-5
RECTANGULAR CARBOY