Datenblatt

Technical data sheet

Sorte | grade: **4b**

Art. Nr. | item No.: **309**

Eigenschaften | characteristics: General purpose filter paper, plain

Anwendung | applications: Suitable for seed testing, germination capacity test and filtration of fine precipitates

Prüfmerkmal property	Prüfmethode test method	Einheit unit	charakteristische Werte typical values
Flächenbezogene Masse Basis weight	DIN EN ISO 536	g/m²	75
Filtriergeschwindigkeit Filtration speed	DIN SPEC 55700	s/10 ml	22,0
Nassberstwiderstand Wet burst strength	DIN ISO 3689	kPa	≥ 30
Wasseraufnahme Water absorption	intern internal	%	270

25.01.2024 Rev.0

All indicated data to be understood as typical average values.
This product is manufactured in conformance with the Quality Management standard DIN EN ISO 9001.

Contact office:

Ahlstrom Germany GmbH
Niederschlag 1
09471 Bärenstein, Germany
T: +49 (0) 37347 / 83-0
F: +49 (0) 37347 / 83-64
filtrak@ahlstrom.com

www.ahlstrom.com

DISCLAIMER: The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as Ahlstrom deems appropriate. Refer to www.ahlstrom.com for contact information.
© Ahlstrom 2023