Micro & Conical Centrifuge Tubes


Available in natural and pure black for light sensitive samples

Microcentrifuge Tubes:
0.2, 0.5 and 1.5 mL *(a-d)*

- Attached flat caps for easy one-handed use
- 0.2 mL is suitable for thermal cyclers
- 0.5 and 1.5 mL versions have a frosted writing area and a syringe pierceable cap
- 100% polypropylene

AUTOCLAVABLE 121° C

NON STERILE



Microcentrifuge Tubes


	Item No.	Color	Vol. *mL*	Outer Dimension *mm*	Closed Height *mm*	Graduations *mL*	Max RCF Rating	Packaging
a	**HS4423**	Natural	0.2	6.9	22.2	None	15,000 x g	1,000/Bulk
b	**HS4422**	Natural	0.5	10.1	31.7	0.1, 0.3, 0.5	15,000 x g	1,000/Bulk
c	**HS4323**	Natural	1.5	12.9	40.2	0.1, 0.5, 1	16,400 x g	500/Bulk
c	**HS4325**	Natural	1.5	12.9	40.2	0.1, 0.5, 1	16,400 x g	1,000/Bulk
d	**HS4323K**	Black	1.5	13.4	40.6	0.1, 0.5, 1, 1.5	18,000 x g	500/Bulk

Conical Centrifuge Tubes:
15 and 50 mL *(e-i)*

- Lot-to-lot traceability with easy tear packaging
- Writing area and flat polyethylene lids for easy labelling
- Imprinted graduations and a self-standing version (only in clear bulk pack)
- Recyclable chipboard used for the racked versions, making them more environmentally friendly

AUTOCLAVABLE 121° C

STERILE R



Skirted

Certified DNase, RNase & Pyrogen Free *(f & h)*

Conical Centrifuge Tubes

	Item No.	Color	Vol. *mL*	Outer Dimension *mm*	Closed Height *mm*	Graduations *mL*	Max RCF Rating	Packaging *500 per Case*
e	**HS4426B**	Natural	15	17	120	0.5	11,000 x g	Bulk/10 x 50pk
	HS4426R	Natural	15	17	120	0.5	11,000 x g	Rack/10 x 50pk
f	**HS4428**	Black	15	17	120	0.5	9,400 x g	Bulk/10 x 50pk
g	**HS4427B**	Natural	50	29	114	5.0	9,500 x g	Bulk/20 x 25pk
	HS4427R	Natural	50	29	114	5.0	9,500 x g	Rack/20 x 25pk
h	**HS4429**	Black	50	29	114	5.0	9,400 x g	Bulk/20 x 25pk
i	**HS4427S**	Natural	50	29 - skirted	114	5.0	9,500 x g	Bulk/20 x 25pk