Dichloromethane

Baker Analyzed

Stabilized with 20 ppm Amylene

Material No: 7053.9025
Effective from: 2013/01/14

# Product Specification

Meets ACS Specifications.

| Test | Specification |
|---|---|
| Assay (by GC) (corrected for water) | >= 99.5 % |
| Appearance | clear |
| Color (APHA) | <= 10 |
| Free Halogens | passes test |
| Preservative (as Amylene) | about 20 ppm |
| Residue after Evaporation | <= 0.002 % |
| Titrable Acid (meq/g) | <= 0.0002 |
| Water (H2O) | <= 0.02 % |
| Trace Impurities (in ppm) | |
| Aluminium (Al) | <= 0.5 |
| Boron (B) | <= 0.02 |
| Barium (Ba) | <= 0.1 |
| Calcium (Ca) | <= 0.5 |
| Cadmium (Cd) | <= 0.05 |
| Cobalt (Co) | <= 0.02 |
| Chromium (Cr) | <= 0.02 |
| Copper (Cu) | <= 0.02 |
| Iron (Fe) | <= 0.1 |
| Magnesium (Mg) | <= 0.1 |
| Manganese (Mn) | <= 0.02 |
| Nickel (Ni) | <= 0.02 |
| Lead (Pb) | <= 0.1 |
| Tin (Sn) | <= 0.1 |
| Zinc (Zn) | <= 0.1 |

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000

Material No: 7053.9025

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000