BRAND. For lab. For life.®

BRAND

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany

BRAND GMBH + CO KG
P.O. Box
97861 Wertheim
Germany
T +49 9342 808 0
info@brand.de
www.brand.de

Product information

Wertheim, 15. Apr. 2024

Product: **Pasteur pipettes PE LD**

Cat. No: **7477 60**

Lot No.: **21510358**

Manufacturing date: **01/2022**

Expiry date: **12/2026**

Very good reproducibility of the number of drops per milliliter. Ideal for aliquots.

The integrated pipetting bulb depresses easily, minimizing fatigue from frequent pipetting.

They are resistant to gas or gamma radiation sterilization.

Graduation [ml]	Withdraw volume with ball [ml]	Outer-Ø tip [mm]	Length [mm]	Drop quantity of drops/ ml	Pack of	Cat. No.
0.5/1.0/1.5/2.0/2.5/3.0	6.2	3	155	20	1000	747760

Best regards,

BRAND GMBH + CO KG

i. A.

Sophie Heller
Product Specialist Repetitive Pipettes and PCR

i. A.

Dr. Nicolas Frank
Product Manager Cell Culture

BRANDGROUP