HIRSCHMANN®

Data sheet

Article: Measuring cylinder 500 ml, low form, class B
Code no.: 240 01 90

Material:
DURAN (Borosilicat glass 3.3)

Chemical properties:
- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Nominal volume:
500 ml

Subdivision:
10 ml

Tolerance:
± 5,0 ml

Height:
max. 255 mm

Standard:
ISO 4788

Hirschmann Laborgeräte GmbH & Co.KG

Eberstadt, 17.12.2021

This letter was created by machine and is valid even without signature.