interscience

BagFilter® Pipet

Lateral filter bag for lab blender with easy pipet access

Ref 111 710 (400 mL)



Standards:

ISO 7218 | ISO 6887-1 | FDA BAM Bacteriological Analytical Manual

Manufactured under:

MADE IN FRANCE

Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagFilter® Pipet is a blender bag with a lateral non-woven filter for all microbiological analyses. A reclosable sticker allows easy pipetting with a micropipette cone. Cross contamination risks are reduced with pipet cones for pathogens pre-enrichment bags. The lower pipetting height reduces risks of MSD (Musculoskeletal risks). It is ideal for homogenizing fibrous samples and for pipetting. The filtration is instantaneous without cross-contamination.

TECHNICAL SPECIFICATIONS

- Bag with lateral filter
- **Reclosable sticker for easy pipetting with a micropipette cone**
- Compatible with **BagRack®** ***Slide*** (Ref 221 440) with BagClip® (Ref 221 040)
- Reduces risks of cross contamination with pipet cones for pathogens pre-enrichment bags
- Reduces risks of Musculoskeletal Disorders (MSD) with a lower pipetting height
- Ideal for homogenizing fibrous samples
- **Multilayer®**: reinforced multicoated complex
- Lateral non-woven filter
- Filter porosity: < 250 microns
- Rigid and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenization of the sample
- Compatible with any lab blender
- International patent pending
- Approved for food contact: Regulation (EC) No 1935/2004
- Gamma ray treated: Gamma 5 to 12 kGy, with certificate
- Average shelf-life: 5 years
- Storage: in a temperate and dry place, away from humidity and light
- Available in 400 mL
- In compliance with ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

BagFilter® Pipet 400

Ref 111 710

- Max blending volume: 400 mL
- Optimal blending volume: 50 - 300 mL
- Optimal pipetting volume : 250 mL
- Bag dimensions: 190 x 300 mm
- Box dimensions: 38 x 22.5 x 15.5 cm, weight: 5.3 kg
- Pack of 25
- Box of 500

BagFilter® Pipet_FT_1118 Pictures and information are not contractually binding. BagFilter® is a trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

www.interscience.com

interscience
30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33 (0)1 34 62 62 61
Email: info@interscience.com

interscience USA & CANADA
32 Cummings Park - Woburn, MA 01801 - USA
P: +1 781 937 0007 - F: +1 781 937 0017
Email: sales.usa@interscience.com

interscience CHINA
上海市徐汇区肇嘉浜路798号坤阳大厦1903室-200030
电话: +86 (0)21-64739390 - +86 189 3097 0733
邮 址: sales.china@interscience.com

interscience SOUTH-EAST ASIA
The Centropod, 80 Changi Rd - 05-07 SINGAPORE 419715
T: +65 6909 0825 - M: +65 8118 5178
Email: sales.asia@interscience.com