
BIOFIL®	Document Number	Revision	Page	Date
	PPT001	A0	1 of 1	10-27-2020
PRODUCT DATA SHEET	Pipette Micro Tips			



Pipette Micro Tips

Description: Pipette tips, natural, Dnase/Rnase-free, Non-pyrogenic, with Filter

Purpose: Use for liquid transfer

Materials


Tips: PP (Polypropylene)
Rack: PP (Polypropylene)

Features

- Preferred accessories for most brand micropipettor
- Tips with PP filter are also available
- Color-coded racks for easy volume identification
- Every rack or case is printed with lot No. for quality traceablity
- DNase/RNase-free and non-pyrogenic
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

② STERILE | R

Dimensions (Unit: mm)



Pipette Micro Tips

Cat. No.	Description	Sterilization	Packaging Configuration
PPT151010	10µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT150010	10µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT451010	10µl, Long, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT450010	10µl, Long, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT151020	20µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT150020	20µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT151100	100µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT150100	100µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT153200	200µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT152200	200µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT451300	300µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT450300	300µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT151000	1000µl, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT150000	1000µl, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case
PPT371000	1000µl, Long, natural, sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	Yes	96 tips/rack, 10 racks/box, 2 boxes/case
PPT370000	1000µl, Long, natural, non-sterilized, Dnase/Rnase-free, Non-pyrogenic, with Filter	No	96 tips/rack, 10 racks/box, 2 boxes/case

Drawn	Approved
Date 27 **Oct** 2020	Date 27 **Oct** 2020
Name Alex Fang	Name James