Ethyl Acetate
Baker Analyzed

Material No: 8037.9025
Effective from: 2017/03/28

Product Specification

Meets ACS Specifications.

Test	Specification
Assay (by GC)	>= 99.5 %
Color (APHA)	<= 10
Residue after Evaporation	<= 0.003 %
Substances Darkened by H2SO4	passes test
Titrable Acid (meq/g)	<= 0.0009
Water (H2O)	<= 0.2 %
Trace Impurities (in ppm)	
Aluminium (Al)	<= 0.500
Barium (Ba)	<= 0.100
Boron (B)	<= 0.020
Cadmium (Cd)	<= 0.050
Calcium (Ca)	<= 0.500
Chromium (Cr)	<= 0.020
Cobalt (Co)	<= 0.020
Copper (Cu)	<= 0.020
Iron (Fe)	<= 0.100
Lead (Pb)	<= 0.100
Magnesium (Mg)	<= 0.100
Manganese (Mn)	<= 0.020
Nickel (Ni)	<= 0.020
Tin (Sn)	<= 0.100
Zinc (Zn)	<= 0.100

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000

Material No: 8037.9025

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000