Valves made of PTFE, two-way and three-way



Two-way / three-way valves made of solid PTFE with hose nozzles. Resistant to chemicals. Very easy connection thanks to its excellent antifriction properties. Blue handles made of PP.

- Resistant to compression max. 1 bar
- Vacuum suitable up to 5 mmHg
- Autoclavable

Valves made of PTFE, two-way and three-way				
Model	For dia. mm	Nominal width mm	Item no.	
three-way valve	4.5	2	8610-0001	N4
three-way valve	6.8	3	8610-0003	N4
three-way valve	10	4	8610-0006	N4
two-way valve	4.5	2	8610-0010	N4
two-way valve	6.8	3	8610-0012	N4
two-way valve	10	4	8610-0015	N4

Hose tubing stop-cock



Simple reflux nozzle with two-sided hose nozzle.

- HDPE natural
- Resistant to acids and chemicals
- Pressure resistant up to 0.3 bar

Hose tubing stop-cock				
For dia. mm	Nominal width mm	Length mm	Item no.	
9-11	7	85	8614-0007	N4
12-15	9	97	8614-0009	N4