IVD In Vitro Diagnostic Medical Device
For professional use only

TECHNICAL DATA SHEET
CASSETTES MACROSTAR V Q PATH

REFERENCES

Cat. No	Description
720-2215	Q Path Macrostrar V white
720-2216	Q Path Macrostrar V green
720-2217	Q Path Macrostrar V pink
720-2218	Q Path Macrostrar V yellow
720-2219	Q Path Macrostrar V blue
720-2220	Q Path Macrostrar V grey
720-2221	Q Path Macrostrar V orange
720-2222	Q Path Macrostrar V red
720-2223	Q Path Macrostrar V lavender



DESCRIPTION

Q Path Macrostrar III with welded lid, large holes.
Its use is reserved for professionals.
To carry and hold operating parts throughout the stages of the histological process.

GENERAL FEATURES

CASSETTE	Measurements (external)	Tolerances
Lenght	83,6 mm	+/- 0,2 mm
Width	28,6 mm	+/- 0,2 mm
Heigth	6 mm	+/- 0,2 mm
Weigth	3,6 gr	+/- 0,2 gr

COMPOSITION

Stain	Polyoxymethylene (POM)
< 5 %	> 95%

CE 05/2010

VWR International SAS
201 rue Carnot
94126 Fontenay-sous-bois
France

IVD In Vitro Diagnostic Medical Device
For professional use only

TECHNICAL DATA SHEET
CASSETTES MACROSTAR V Q PATH

PACKAGING

Type of Packaging	Dimensions	Quantity	Total Quantity
Inner Box	240 x 170 x 310	3 x 500	1500
Outer Box	531 x 244 x 327	1500	1500
Pallet	1200 x 800 x 1890	30 x 1500	45000

APPLICABLE STANDARDS

- Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices

CE 05/2010

VWR International SAS
201 rue Carnot
94126 Fontenay-sous-bois
France