9F, No. 108, Jhouzih St., Taipei, Taiwan
Tel: + 886-2-8751-1888 Fax: + 886-2-6602-1218
E-mail: sales@abnova.com

Datasheet

IgG (Mouse) Isotype Control (LE/AF)

Catalog Number: PAB19464

Regulation Status: For research use only (RUO)

Product Description: Mouse IgG raised from normal mouse serum.

Host: Mouse

Applications: ELISA, Flow Cyt
(See our web site product page for detailed applications information)

Protocols: See our web site at http://www.abnova.com/support/protocols.asp or product page for detailed protocols

Form: Liquid

Conjugation: LE/AF

Concentration: 0.5 mg/mL

Isotype: IgG

Recommend Usage: Flow Cytometry (? 1 ug/10^6 cells)
ELISA (1:500-1:5000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS

Storage Instruction: Store at -20°C. For long term storage store at -80°C.
Aliquot to avoid repeated freezing and thawing.